





**FX.1914.REF FIXADOR CLICK FIX MEDIO  
(PARA MEMBROS SUPERIORES)**



- Produto Estéril
- Para Pinos de  $\varnothing 4.0$ ,  $\varnothing 3.5$ ,  $\varnothing 3.0$
- Barras Longitudinal Alumínio  $\varnothing 8.0$  com comprimento de 100, 125, 150, 175, 200, 250 e 300mm
- Barras Transversal Inox  $\varnothing 4,0$  com comprimento de 100, 125, 150 e 200mm



**ENCAIXE RÁPIDO POR FRENTE**